

ELECTRODE ACCESSORIES


Accessories

Jumper Cables

The electrode jumper cables allow the same electrode to be plugged into two different inputs on the preamp. Individually packaged. (1.5mm connector)

Red — 301003

Black — 301004



Electrode Extension Cables

The electrode extension cable plugs into any 1.5mm connector electrode adding 2 meters to any existing lead. Sold individually – multi-colored cables available.

301005



Electrode Adapters

The electrode adapters are used to convert between the 2mm pin-style connector and the 1.5mm touchproof safety connector.

Description	Part No.
2mm Female to 1.5mm Female ("old style" 2mm pin electrodes to new 1.5mm preamp)	301201
1.5mm Male to 2mm Male (new 1.5mm electrodes to old 2mm preamp) RED Adapter	301202
1.5mm Male to 2mm Male (new 1.5mm electrodes to old 2mm preamp) BLACK Adapter	301203
1.5mm Male to 2mm Male (new 1.5mm electrodes to old 2mm preamp) GREEN Adapter	301204
1.5mm Male to 2mm Male (new 1.5mm electrodes to old 2mm preamp)	301205



1.5mm Male to 2mm Male



2mm Female to 1.5mm Female

To Order Call 1-800-456-5205

Bio-logic
SYSTEMS CORP