


FREQUENTLY USED PRODUCTS

Disposable Electrodes

Disc with Lead


These pre-gelled, silver/silver chloride electrodes provide high quality recordings with the convenience of a disposable. Bio-logic electrodes are ideal for SEP, AEP and ENG studies. Electrodes with a lead cable end with a 1.5mm connector. *SINGLE PATIENT USE ONLY.*



Length	Part No.
.75m — Lead Color: Black, Red, Green, White and Yellow 1 Card (10/pkg)	101623
1.5m — Lead Color: Black, Red, Green and Blue 15 Cards (60/pkg)	101621
2.5m — Lead Color: Black, Red, Green and Blue 15 Cards (60/pkg)	101622

Ground with Lead

The pre-gelled silver/silver chloride ground electrodes have a large recording area with a black or green attached lead (12 green and 12 black/pkg). The electrodes are available in two different lengths: 2 meter and 2.5 meter. These electrodes are great for IOM and Sleep recordings. *SINGLE PATIENT USE ONLY.*



Lead Length	Part No.
2m (24/pkg)	101620
2.5m (24/pkg)	101625

Skin Preparation Material

NuPrep Skin Preparation

Topical skin abrasive. For low impedance. Mild abrasive for patient comfort. Chloride and acetone free.

109090



DVD/PD-CD

Optical Disks Sold Individually

Description	Part No.
4.7 GB DVD (single-side) Optical Disk	121247
5.2 GB DVD Optical Disk	121252
9.4 GB DVD Optical Disk	121294
650 MB PD/CD Optical Disk	121265



Gels and Pastes

Ten20 Conductive Paste

Non-drying, water-soluble conductive paste. Recommended for EP and EEG disc electrodes and also used under collodion electrodes.

Description	Part No.
4 oz. Jar	109100
4 oz. Tube	109101
8 oz. Jar	109102



To Order Call 1-800-456-5205

Bio-logic
SYSTEMS CORP