


SURFACE ELECTRODES - DISC

Disc Electrodes

Reusable

Bio-logic offers the widest variety of reusable disc electrodes, all ending with a molded connector. Adult and pediatric size disc electrodes are available in the following cable lengths: 1m, 1.5m and 3m (adult only). The silicone rubber leads (latex free) are flexible and tangle-free. The cast discs are manufactured from "fine" 99.99% pure silver. All electrodes end with a 1.5mm connector. Multi-colored leads. (10/pkg)



Adult - 9mm Cast Disc

Silver Disc	Part No.
1m Silicone Cable	101041
1.5m Silicone Cable	101042
3m Silicone Cable	101043

Silver Chloride Disc	Part No.
1m Silicone Cable	101051
1.5m Silicone Cable	101052
3m Silicone Cable	101053

Gold Disc	Part No.
1m Silicone Cable	101061
1.5m Silicone Cable	101062
3m Silicone Cable	101063

Tin Disc	Part No.
1m Teflon Cable	101071
1.5m Teflon Cable	101072
3m Teflon Cable	101073

Pediatric - 6mm Cast Disc

Silver Disc	Part No.
1m Silicone Cable	101044
1.5m Silicone Cable	101045

Silver Chloride Disc	Part No.
1m Silicone Cable	101054
1.5m Silicone Cable	101055

Gold Disc	Part No.
1m Silicone Cable	101064
1.5m Silicone Cable	101065

Tin Disc	Part No.
1m Teflon Cable	101074
1.5m Teflon Cable	101075

All Leads in 1.5mm Connector

To Order Call 1-800-456-5205

Bio-logic[®]
SYSTEMS CORP

SURFACE ELECTRODES – DISC

Disc Electrodes – Reusable

Overmolded 9mm Disc


Lightweight electrode features a high quality stamped disc, which is overmolded with a plastic cover. The plastic cover provides excellent strain relief for a long electrode life and an insulator to prevent loss of gel/paste. Ten multi-colored 1 meter leads ending in a 1.5mm molded connector.

1m Cable	Part No.
Silver	101081
Gold	101087



High Dome 10mm Disc

The High Dome (Grass-style) electrodes are available in two different metals: gold and silver as well as two styles of cables: silicone (tangle-free) and Teflon (light-weight). Ten multi-colored leads ending in a 1.5mm molded connector.



Description	Part No.
1m Silver Teflon Cable	100001
1m Gold Teflon Cable	100021
1m Silver Silicone Cable	100041
1m Gold Silicone Cable	100061
1.5m Gold Teflon Cable	100022
1.5m Gold Silicone Cable	100062

Disposable EEG/EP Cup Electrode

Two silver/silver chloride disposable EEG/EP cup electrodes are now available. Simply snap the electrode into the reusable cable and proceed as usual. There is no messy clean-up. After testing, remove the disc from the cable and discard. Like the traditional cup electrodes, the disposable electrodes can be used with conductive paste or collodion. The reusable lead cables are 1 meter in length and end in a 1.5mm connector.



Standard



Bristle

Description	Part No.
Standard (30 qty)	101653
Standard (500 qty)	101655
Bristle (30 qty)	101654
Bristle (500 qty)	101656
Reusable Lead Cables 1m (10 qty)	301650



Reusable Lead Cable

To Order Call 1-800-456-5205


Bio-logic
SYSTEMS CORP

SURFACE ELECTRODES

Disposable Electrodes

Disc with Lead


These pre-gelled, silver/silver chloride electrodes provide high quality recordings with the convenience of a disposable. Bio-logic electrodes are ideal for SEP, AEP and ENG studies. Electrodes with a lead cable end with a 1.5mm connector. *SINGLE PATIENT USE ONLY.*



Length	Part No.
.75m — Lead Color: Black, Red, Green, White and Yellow 1 Card (10/pkg)	101623
1.5m — Lead Color: Black, Red, Green and Blue 15 Cards (60/pkg)	101621
2.5m — Lead Color: Black, Red, Green and Blue 15 Cards (60/pkg)	101622

Snap Electrodes

Custom-designed electrodes with soft cloth and rounded corners for patient comfort. They use a 40mil hypoallergenic hydrogel for better adhesion and electrical performance yet are repositionable. They provide optimum performance and comfort to meet the special needs of your newborn and adult patients. Three-pack foil pouches keep them fresh. *SINGLE PATIENT USE ONLY.*



Shape	Part No.
Rectangle (60/pkg)	101600
Rectangle (300/pkg)	101603

Rectangular Snap Electrode
(1.31" x 8.75") (3.5 x 2.4cm)

Snap Electrodes

The square snap electrode uses a specially formulated adhesive solid gel with 100 percent adhesive surface area. Use this electrode when aggressive adhesion is desirable. *FOR SINGLE PATIENT USE ONLY. NOT RECOMMENDED FOR USE ON INFANTS.*



Shape	Part No.
Square (30/pkg)	101604
Square (300/pkg)	101605

Square Snap Electrode
(1" x 1") (2.6 x 2.6cm)


To Order Call 1-800-456-5205

SURFACE ELECTRODES

Disposable Electrodes

Ground with Lead

The pre-gelled silver/silver chloride ground electrodes have a large recording area with a black or green attached lead (12 green and 12 black/pkg). The electrodes are available in two different lengths: 2 meter and 2.5 meter. These electrodes are great for IOM and Sleep recordings. *SINGLE PATIENT USE ONLY.*



Lead Length	Part No.
2m (24/pkg)	101620
2.5m (24/pkg)	101625

Clip-Type

Pre-gelled silver/silver chloride electrodes provide a high quality recording with the convenience of a disposable. Made with a conductive hydrogel, these electrodes can be reapplied multiple times to the same patient. Packaged 10 electrodes per card in a shrink wrap to prevent drying out. For use with reusable lead clips. (100/pkg) *SINGLE PATIENT USE ONLY.*

Description	Part No.
Ground Electrode (100/pkg)	101632
Ring Electrode (100/pkg)	101633
Disc Electrode (100/pkg)	101634
Tab Electrode (100/pkg)	101635



101634
.5" square (1.5cm)



101635
.875" x 1"
(2.2cm x 2.5cm)



101632
1.75" x 1.25"
(4.5cm x 3.1cm)



101633
4" x .25" (10cm x .6cm)

To Order Call 1-800-456-5205


Bio-logic[®]
SYSTEMS CORP

SURFACE ELECTRODES

Disposable Electrode Lead Cables

Tabloc

The low profile Tabloc lead easily slides on the tab of a disposable surface electrode. Each package contains 5 different colored clips (red, green, black, white and blue) and is available in two different lengths.



Length	Part No.
1m	301645
3m	301646

Pinch Clip


The Pinch Clip lead easily attaches to a disposable snap electrode. 18" lead wires.



Description	Part No.
3 Clips (red, white and blue)	301663
4 Clips (red, white, blue and black)	301664
5 Clips (red, white, blue, green, and black)	301665
6 Clips (red, white, blue, green, black and yellow)	301666

Alligator Clip

The reusable Alligator Clip lead is for use with clip-type disposable surface electrodes. Each package contains 3 (red, green and black) leads.



Length	Part No.
1m	301636
30cm	301638

Snap Lead

For use with disposable snap electrodes. The color-coded snaps (red, green, blue, white and black) are sold in a package of five.



Length	Part No.
1m	301601
1.5m	301604
3m	301602

To Order Call 1-800-456-5205

Bio-logic
SYSTEMS CORP


SURFACE ELECTRODES – EEG, EP & EMG

Ear Clip Electrodes

Reusable

Two disc-type electrodes mounted on a spring clip with an attached 1 meter lead (1.5mm connector). Available in three metal types: Silver, Silver Chloride and Gold and two disc sizes 6mm (pediatric) and 10mm (adult). One red and one blue electrode per package.

Description	Part No.
Silver – Adult	101204
Silver – Pediatric	101203
Gold – Adult	101224
Gold – Pediatric	101223
Silver Chloride – Adult	101214
Silver Chloride – Pediatric	101213



Adult




Pediatric

Disposable

Lightweight adult ear clip electrodes feature detachable silver chloride discs and 1 meter attached lead (1.5mm connector). One red and one blue electrode per package.



Silver Chloride



Replacement Discs

Description	Part No.
Silver Chloride	101234
Replacement Discs (4)	101244

To Order Call 1-800-456-5205

Bio-logic[®]
SYSTEMS CORP

SURFACE ELECTRODES – EEG, EP & EMG

Digital Ring Electrodes

Used for sensory nerve studies, the Ring Electrodes adjust to any size finger with a slide adjustment. Come as a pair with 1 meter attached leads.



Description	Part No.
Stainless Steel	101953
Gold	101923

EMG Disc Electrodes


The EMG Disc Electrodes are round (10mm), flat and made of stainless steel with locating pins. The discs are sold in pairs with untwisted red and black 1 meter leads.

101833



Bar Electrodes

The bar electrode consists of two convex discs mounted 30mm apart and embedded in a plastic bar with an attached 1 meter. The Bipolar Felt Pad electrode uses two felt pads mounted 30mm apart and is ideal for stimulation. Felt tips should be soaked in saline prior to use.



Description	Part No.
Bar Electrode	101803
Bipolar Bar Electrode	101813

Replacement

Felt Pads

(10/pkg) 101899



Removable

T-Handle

101898



Pediatric Stimulating Electrode

The Pediatric Stimulating Electrode has a "double banana" type plug spaced 12.7mm apart and is useful for children and blink reflex testing.

101822



To Order Call 1-800-456-5205

Bio-logic
SYSTEMS CORP


SURFACE ELECTRODES – EEG, EP & EMG

Saline and E.S.D. Grounding Straps

Saline

The wrap around saline ground straps are 20mm wide with a 1.25 meter cable. Available in two sizes, requires soaking in saline.

Description	Part No.
Adult (55cm)	101854
Pediatric (25cm)	101853



Pediatric



Adult

E.S.D.

The E.S.D. Ground straps use conductive Velcro and are available in four sizes. They can be used dry and include a 1 meter lead cable.

Size	Part No.
Small (18cm)	101871
Medium (36cm)	101873
Large (76cm)	101875
X-Large (117cm)	101877



Ground Electrode

The Ground Plate Electrode is a stainless steel round disc (30mm) with a 1 meter attached lead.

101863



To Order Call 1-800-456-5205

Bio-logic[®]
SYSTEMS CORP